

In-solution digestion

<Notice!> Wearing cleaned gloves is needed, and manipulating in laminar flow cabinet is preferred.

Prepare reagents before use:

- ✓ 50 mM ammonium bicarbonate (ABC) solution (pH 8.5)
- ✓ 12.5 ng/ μ L trypsin stock

Procedure

1. Add 50 μ L of 50 mM ABC solution to redissolve proteins from acetone precipitation.
2. Add 48 mg of urea into sample. Final concentration of urea is 8M and the volume is ~100 μ L.
3. Add 1 μ L of 1 M DTT. Incubate at RT for 30 min.
4. Add 5 μ L of 1 M IAM. Incubate at RT for 30 min.
5. Add 400 μ L of 50 mM ABC
6. Add trypsin (w/w, enzyme: protein=1: 20~50) and incubate at 37°C for at least 16 hr.
7. Perform stage tip desalting.
8. Vacuum dry the sample.
9. Add 0.1% formic acid.